

Product Datasheet

Anti-CD20/PE-Cy5.5 (orb971935)

Description	Anti-CD20/PE-Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD20
Target	human CD20
Preservatives	0.01M PBS, 0.5% rAlbumin, 0.02% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Uniprot ID	P11836
Expiration Date	12 months from date of receipt.