

Product Datasheet

**Anti-Phospho-TrkA(Tyr674 + Tyr675) +
TrkB(Tyr706 + Tyr707)/AP (orb972860)**

Description

Anti-Phospho-TrkA(Tyr674 +

Tyr707)/AP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

AP

Tested

ELISA, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human TrkA around the phosphorylation site of Tyr674/675

Target

Phospho-TrkA(Tyr674 + Tyr675) + TrkB(Tyr706 + Tyr707)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

P04629

Dilution Range

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.