



Product Datasheet

PGC-1alpha Antibody (orb97371)


Description	Mouse monoclonal antibody to PGC-1alpha
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human PGC-1alpha expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	91kDa
Entrez	10891
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 (Lane1) ...



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human PGC-1alpha...