

Product Datasheet

S6 Ribosomal Protein (phospho-Ser235+Ser236) antibody (FITC) (orb9738)

Description

Rabbit polyclonal to S6 Ribosomal

antibody conjugated to FITC

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

FITC

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human S6 Ribosomal Protein around the phosphorylation site of Ser235/236

Target

RPS6

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

27kDa

Uniprot ID
P62753
Dilution Range

IF: 1:50-200

Expiration Date

12 months from date of receipt.