



Product Datasheet

MSH6 Antibody (orb97388)


Description	Mouse monoclonal antibody to MSH6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of human MSH6 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	160kDa
Entrez	2956
Dilution Range	WB: 1/500 - 1/2000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human MSH6 (AA:...



Line graph illustrates about the Ag-Ab r...



Western blot analysis of MCF-7 (Lane1), ...