


Product Datasheet


RSK2 Antibody (orb97391)

Description	Mouse monoclonal antibody to RSK2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Purified recombinant fragment of human RSK2 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	80kDa
Entrez	6197
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.




Flow
cytometric
analysis of
HepG2 cells

...



Western blot
analysis of
Hela (Lane1),
M...



Western blot
analysis of
human RSK2
(AA:...