



## Product Datasheet

### GSTM1 Antibody (orb97400)


|                            |                                                                             |
|----------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to GSTM1                                          |
| <b>Species/Host</b>        | Mouse                                                                       |
| <b>Reactivity</b>          | Human, Monkey                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                |
| <b>Tested Applications</b> | ELISA, FC, IHC, WB                                                          |
| <b>Immunogen</b>           | Purified recombinant fragment of human GSTM1 expressed in E. Coli.          |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                                |
| <b>Concentration</b>       | 1mg/ml                                                                      |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.                 |
| <b>Note</b>                | For research use only                                                       |
| <b>Isotype</b>             | Mouse IgG1                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                  |
| <b>MW</b>                  | 26kDa                                                                       |
| <b>Entrez</b>              | <a href="#">2944</a>                                                        |
| <b>Dilution Range</b>      | WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                             |



Flow  
cytometric  
analysis of  
Jurkat cells...



Western blot  
analysis of  
Cos7 (Lane1),  
M...



IHC-P of liver  
cancer  
tissues using  
GSTM...