



Product Datasheet

CD9 Antibody (orb97403)


Description	Mouse monoclonal antibody to CD9
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Synthesized peptide of human CD9.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	25kDa
Entrez	928
Dilution Range	IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



IHC-P of kidney cancer tissues using CD9...



Western blot analysis of HEK293 (Lane1) ...



Line graph illustrates about the Ag-Ab r...