



Product Datasheet

SOX9 Antibody (orb97410)


Description	Mouse monoclonal antibody to SOX9
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Purified recombinant fragment of human SOX9 expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	56kDa
Entrez	6662
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human SOX9 (AA:...)



Western blot analysis of Lovo (Lane1) and another sample (Lane2)



Confocal immunofluorescence analysis of ...