



Product Datasheet

SERPINE1 Antibody (orb97411)


Description	Mouse monoclonal antibody to SERPINE1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Purified recombinant fragment of human SERPINE1 expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	45kDa
Entrez	5054
Dilution Range	IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
NIH/3T3
cell...



IHC-P of
kidney
cancer
tissues using
SER...



Western blot
analysis of
HEK293
(Lane1) ...