



Product Datasheet

troponin T2 Antibody (orb97466)


Description	Mouse monoclonal antibody to troponin T2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human troponin T2 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	36kDa
Entrez	7139
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human (AA: 88-2...



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of MCF-7 cells

...