

Product Datasheet

Anti-phospho-Smad2/Smad3 (Thr8)/Cy3 (orb974699)

Description	Anti-phospho-Smad2/Smad3 (Thr8)/Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Cy3
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Smad2/Smad3 around the phosphorylation site of Thr8
Target	phospho-Smad2/Smad3 (Thr8)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15796
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.