



Product Datasheet

DLL4 Antibody (orb97478)


Description	Mouse monoclonal antibody to DLL4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human DLL4 expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b/4A11F8 Mouse IgG1
Clonality	Monoclonal
MW	74.6kDa
Entrez	54567
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HeLa cells u...



Line graph
illustrates
about the Ag-
Ab r...



Western blot
analysis
using DLL4
HEK293 ...