



Product Datasheet

LPA Antibody (orb97485)


Description	Mouse monoclonal antibody to LPA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	Purified recombinant fragment of human LPA (AA: 4330-4521) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	501kDa
Entrez	4018
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



IHC-P of breast cancer tissues using LPA...



Confocal immunofluorescence analysis of ...



Western blot analysis of human recombina...