



Product Datasheet

E2F1 Antibody (orb97489)


Description	Mouse monoclonal antibody to E2F1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Purified recombinant fragment of human E2F1 (AA: 69-223) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	46.9kDa
Entrez	1869
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



IHC-P of rectum cancer tissues using E2F...



Western blot analysis of human recombina...



Flow cytometric analysis of HeLa cells u...