



Product Datasheet

ARHGDI2 Antibody (orb97492)


Description	Mouse monoclonal antibody to ARHGDI2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Purified recombinant fragment of human ARHGDI2 (AA: FULL(1-204)) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	26kDa
Entrez	396
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western blot analysis of human recombinant...



Western blot analysis of Jurkat (Lane1),...