


## Product Datasheet

### IKBKB Antibody (orb97514)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IKBKB                           |
| <b>Species/Host</b>        | Mouse                                                        |
| <b>Reactivity</b>          | Human                                                        |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | ELISA                                                        |
| <b>Immunogen</b>           | Purified recombinant fragment of IKBKB expressed in E. Coli. |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                 |
| <b>Concentration</b>       | 1mg/ml                                                       |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.  |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | Mouse IgG1                                                   |
| <b>Clonality</b>           | Monoclonal                                                   |
| <b>MW</b>                  | 43kDa                                                        |
| <b>Entrez</b>              | <a href="#">3551</a>                                         |
| <b>Dilution Range</b>      | ELISA: 1/10000                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |



Western blot analysis of recombinant pro...