


## Product Datasheet

### IKBKE Antibody (orb97534)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IKBKE                                     |
| <b>Species/Host</b>        | Mouse                                                                  |
| <b>Reactivity</b>          | Human                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | ELISA                                                                  |
| <b>Immunogen</b>           | Purified recombinant fragment of IKBKE (aa1-257) expressed in E. Coli. |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                           |
| <b>Concentration</b>       | 1mg/ml                                                                 |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.            |
| <b>Note</b>                | For research use only                                                  |
| <b>Isotype</b>             | Mouse IgG1                                                             |
| <b>Clonality</b>           | Monoclonal                                                             |
| <b>Entrez</b>              | <b>9641</b>                                                            |
| <b>Dilution Range</b>      | ELISA: 1/10000                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Western blot analysis of recombinant pro...