



Product Datasheet

Glucose-6-phosphate isomerase Antibody (orb97545)


Description	Mouse monoclonal antibody to Glucose-6-phosphate isomerase
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	Purified recombinant fragment of human GPI expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	63kDa
Entrez	2821
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells t...



Confocal immunofluorescence analysis of ...



Immunohistochemical analysis of formalin...