



Product Datasheet

PYK2 Antibody (orb97549)

Description	Mouse monoclonal antibody to PYK2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Purified recombinant fragment of PYK2 (aa815-997) expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	115.8kDa
Entrez	2185
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of Raji (Lane1), P...