



Product Datasheet

PEG10 Antibody (orb97563)


Description	Mouse monoclonal antibody to PEG10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Purified recombinant fragment of human PEG10 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Isoform RF1 (37kDa); RF1/RF2 (80kDa)
Entrez	23089
Dilution Range	IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of HepG2 (Lane1), ...



IHC-P of human Placenta tissues using PE...