



Product Datasheet

Calnexin Antibody (orb97585)

Description	Mouse monoclonal antibody to Calnexin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Synthetic peptide corresponding to aa (CEAAEERPWLWVYILTVAL) of human Calnexin, conjugated to KLH.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	90kDa
Entrez	821
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 (Lane1), H...



Confocal immunofluorescence analysis of ...