


Product Datasheet

HSPA13 antibody (orb97642)

Description	Goat polyclonal antibody to HSPA13
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HSPA13 / STCH (aa83-95)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in Human Uterus lysates (calculated MW of 51.9kDa according to NP_008879.3).The protein is commonly known as a 60 or 70kDa heat shock protein (see nomenclature).Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	ESVELADSNPQNT
Clonality	Polyclonal
MW	51.9
Entrez	6782
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Uterus ly...