



Product Datasheet

Stomatin antibody (orb97643)

Description	Goat polyclonal to Stomatin
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Stomatin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Human Erythrocytes, Spleen, Bone Marrow and Lung lysates (calculated MW of 31.7kDa according to NP_004090.4). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	EKRHTRDSEAQR
Clonality	Polyclonal
MW	31.7
Entrez	2040
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Li...



Western blot analysis of Human Lung lysis...