




Product Datasheet

MET antibody (orb97647)

Description	Mouse monoclonal antibody to MET
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	MET
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	155541
Uniprot ID	P08581
NCBI	NP_001120972.1, NP_000236.2
Dilution Range	IF/ICC: 1:100, WB: 1:100-500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of HepG2 cell line...