




Product Datasheet

LC3 antibody (orb97657)

Description	Mouse monoclonal antibody to LC3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	This LC3 antibody is generated from mouse immunized with a full length recombinant protein of human LC3 (APG8).
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1 k
Clonality	Monoclonal
Uniprot ID	Q9H492 , Q9GZQ8
NCBI	NP_115903.1 , NP_852610.1
Dilution Range	WB: 1:1000, IF/ICC: 1:200, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cell lysat...



Immunofluorescence analysis of U251 cell...