




Product Datasheet

ALDH1A1 antibody (orb97664)

Description	Mouse monoclonal antibody to ALDH1A1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	ALDH1A1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	54862
Uniprot ID	P00352
NCBI	NP_000680.2
Dilution Range	WB: 1:500-1000, IF/ICC: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Immunofluorescence analysis of NCI-H460 ...