




Product Datasheet

AKT2 antibody (orb97665)

Description	Mouse monoclonal antibody to AKT2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	AKT2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	55769
Uniprot ID	P31751
NCBI	NP_001617.1
Dilution Range	WB: 1:100-500, IHC-P: 1:50-100, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of mouse brain tis...