

Product Datasheet

HSPA1A antibody (orb97670)

Description	Mouse monoclonal to HSPA1A. This HSPA1A monoclonal antibody is generated from mouse immunized with HS...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	This HSPA1A antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 574-600 amino acids from human HSPA1A.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
Uniprot ID	P08107
NCBI	NP_005336.3, NP_005337.2
Dilution Range	WB: 1:100-2000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.

