




Product Datasheet

c-Kit antibody (orb97685)

Description	Mouse monoclonal antibody to c-Kit
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	KIT
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	109865
Uniprot ID	P10721
NCBI	NP_000213.1 , NP_001087241.1
Dilution Range	IF/ICC: 1:10-50, WB: 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SY5Y cell...



Western blot analysis of A549 cell line ...