


Product Datasheet

hMSRB1 antibody (orb97723)


Description	Anti-human-MSRB1 (h-Methionine sulfoxide reductase B1) is a rabbit polyclonal antibody that was prod...
Species/Host	Rabbit
Conjugation	Unconjugated
Form/Appearance	liquid (Supplied as serum, preserved in glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Anti-human MsrB1 can be used for Western blot, ELISA and Immunocytochemistry.
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using hMS...



Immunofluorescence analysis of using hMS...



Immunofluorescence analysis of using hMS...