

Product Datasheet

Anti-Anti-phospho-Smad2 (Ser250)/Cy3 (orb977339)

Description

Anti-Anti-phospho-Smad2 (Ser250)/Cy3

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Cy3

Tested

FC

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Anti-phospho-Smad2 (Ser250)

Target

SMAD2

Preservatives

Preservative: 15mM Sodium Azide,
Constituents: 1% rAlbumin, 0.01M PBS,
pH 7.4

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

58kDa

Uniprot ID
Q15796
Expiration Date

12 months from date of receipt.