

## Product Datasheet

# Anti-Phospho-Insulin Receptor Beta (Tyr1361)/Cy3 (orb977539)

|                            |                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-Phospho-Insulin Receptor Beta (Tyr1361)/Cy3                                                                              |
| <b>Species/Host</b>        | Rabbit                                                                                                                        |
| <b>Reactivity</b>          | Human                                                                                                                         |
| <b>Conjugation</b>         | Cy3                                                                                                                           |
| <b>Tested Applications</b> | FC                                                                                                                            |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1361 |
| <b>Target</b>              | Phospho-Insulin Receptor Beta (Tyr1361)                                                                                       |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide,<br>Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4                                              |
| <b>Form/Appearance</b>     | Liquid                                                                                                                        |
| <b>Concentration</b>       | 1mg/ml                                                                                                                        |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                              |
| <b>Note</b>                | For research use only                                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                                    |
| <b>CAS Number</b>          | 3                                                                                                                             |
| <b>Uniprot ID</b>          | <b>P06213</b>                                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                               |