

Product Datasheet

SARS-ACN/1 protein (orb97844)

Description	The protein contains the Nucleocapsid protein immunodominant fragments, amino acids: 340-390. SARS-A...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl, 60 mM NaCl and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Recombinant SARS-ACN Antigen may be used in ELISA and Western blots, excellent for detection of SARS with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
MW	Theoretical MW: 32 kDa
Expiration Date	6 months from date of receipt.