

Product Datasheet

HAV-P2C-P3B protein (orb97867)

Description	TheE.coli derived 53.7 kDa recombinant protein contains the P2C-P3B immunodominant regions, amino ac...
Form/Appearance	liquid (Supplied in 10 mM Tris-HCl pH 9.6, 0.1% SDS and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Recombinant HAV-P2C-P3B may be used in ELISA and Western blots, excellent for detection of HAV with minimal specificity problems.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
MW	Theoretical MW: 53.7 kDa
Expiration Date	6 months from date of receipt.