

Product Datasheet

HBV-Core 18 kDa protein (orb97876)

Description	Hepatitis B Virus Core Antigen, molecular weight: 18 kDa, contains the fragments of the Core protein...
Form/Appearance	liquid (Supplied in 1x PBS, 75 mM NaCl and 20% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Immunochromatography (capture and conjugate) Immunogen for monoclonal antibody production ELISA.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, RP-HPLC)
Expiration Date	6 months from date of receipt.