

Product Datasheet

HCV-N-G1b protein (orb97888)

Description	The protein contains the HCV core nucleocapsid immunodominant regions, amino acids: 2-119. Hepatitis...
Form/Appearance	liquid (Supplied in 25 mM Tris-HCl pH 8.0, 1.5 M urea, 0.2% Triton-X and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
MW	Theoretical MW: 40 kDa
Expiration Date	6 months from date of receipt.