

Product Datasheet

HCV-N-G2a protein (orb97889)

Description	The protein contains the HCV core nucleocapsid immunodominant regions, amino acids: 1-102. The prote...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl pH 8.0 and 5 mM EDTA)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
Expiration Date	6 months from date of receipt.