

Product Datasheet

Anti-phospho-TrkA (Tyr490)/Cy3 (orb979216)

Description	Anti-phospho-TrkA (Tyr490)/Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy3
Tested Applications	FC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human TrkA isoform 1(790aa) around the phosphorylation site of Tyr490
Target	phospho-TrkA (Tyr490)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04629
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.