

Product Datasheet

SHP2 (phospho-Ser576) antibody (FITC) [Out of stock] (orb9793)

Description

Rabbit polyclonal to SHP2 (phospho-

Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SHP1 around the phosphorylation site of Ser576
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	139kDa
Uniprot ID	O15357, Q9WVR3, Q6P549
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.