

## Product Datasheet

**Anti-Phospho-Dab1 (Tyr198) /Cy3 [Out of stock]  
(orb980187)**

**Description**

Anti-Phospho-Dab1 (Tyr198) /Cy3

**Species/Host**

Rabbit

**Conjugation**

Cy3

**Tested Applications**

ICC, IF

**Immunogen**

KLH conjugated synthesised phosphopeptide derived from human Dab1 around the phosphorylation site of Tyr198 [V(p-Y)QY]

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

1mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

65kDa

**Uniprot ID**

**O75553**, **P97318**, **Q8CJH2**

**Dilution Range**

ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

**Expiration Date**

12 months from date of receipt.