

Product Datasheet

Rabbit anti Mouse IgG1 IgG2a IgG2b IgG3 (TRITC) antibody (orb98171)

Description	Rabbit polyclonal to IgG1 IgG2a IgG2b IgG3 (TRITC)
Species/Host	Rabbit
Conjugation	TRITC
Storage	Storage: Lyophilized – can be stored at +4°C. After reconstitution aliquots can be stored at -2°C
Note	For research use only
Clonality	Polyclonal
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.
Expiration Date	12 months from date of receipt.