


Product Datasheet

Sheep anti Human IgG1 IgG2 IgG3 IgG4 antibody [Out of stock] (orb98174)

Description	Sheep polyclonal to IgG1 IgG2 IgG3 IgG4
Species/Host	Sheep
Preservatives	Sheep Serum anti Human IgG1 IgG2 IgG3 IgG4 (subclass specific)
Storage	Lyophilized – can be stored at +4°C. After reconstitution aliquots can be stored at -2°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Sheep anti
Human
IgG1 IgG2
IgG3 IgG4
ant...