

Product Datasheet

Agarose Nr.4 (orb98187)

Description	Polyclonal antibody to agarose Nr4 which is an especially formulated agarose to prepare gel blocks f...
Storage	Storage: Product should be stored at 4-8°C. DO NOT FREEZE
Note	For research use only
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.
Expiration Date	12 months from date of receipt.