

Product Datasheet

Anti-phospho-MAP4K1(Ser171)/Cy3 (orb982175)

Description	Anti-phospho-MAP4K1(Ser171)/Cy3
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Cy3
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAP4K1 around the phosphorylation site of Ser171
Target	phospho-MAP4K1 (Ser171)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	91kDa
Uniprot ID	Q92918
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500
Expiration Date	12 months from date of receipt.