

Product Datasheet

Anti-Phospho-Histone H2A.X (Tyr143)/Cy3 (orb982256)

Description	Anti-Phospho-Histone H2A.X (Tyr143)/Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy3
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Histone H2AX around the phosphorylation site of Tyr143
Target	Phospho-H2AX (Tyr143)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16104
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.