

Product Datasheet

Basal cell Cytokeratin antibody (orb98233)

Description

Mouse monoclonal to Basal cell

Species/Host	Mouse
Reactivity	Aves, Canine, Gallus, Guinea pig, Hamster, Human, Porcine, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, WB
Preservatives	Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.
Form/Appearance	Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.
Storage	Storage: The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Source	RCK103 is a Mouse monoclonal IgG1 antibody derived by fusion of SP2/0-Ag14 Mouse myeloma cells with spleen cells from a BALB/c Mouse immunized with a mix of cell preparations containing Human Cytokeratins.
Uniprot ID	multiple