

Product Datasheet

Anti-Phospho-RPS6KA1 (Thr573)/Cy3 [Out of stock] (orb982720)

Description	Anti-Phospho-RPS6KA1 (Thr573)/Cy3 [Out of stock]
Species/Host	Rabbit
Conjugation	Cy3
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6KA1 around the phosphorylation site of Thr573
Target	Phospho-RPS6KA1 (Thr573)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15418
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, IF: 1:100-500, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.