

Product Datasheet

**Anti-phospho-PI3 kinase p85 alpha + gamma
(Tyr467 + Tyr199)/Cy3 (orb982735)**

Description

Anti-phospho-PI3 kinase p85 alpha

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy3
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467
Target	phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27986
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500
Expiration Date	12 months from date of receipt.