

Product Datasheet

Anti-phospho-Bcl-xL (Ser62)/Cy3 (orb983351)

Description	Anti-phospho-Bcl-xL (Ser62)/Cy3
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Cy3
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Bcl-xL around the phosphorylation site of Ser62
Target	phospho-Bcl-xL (Ser62)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q07817
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.